

MULTI PARAMETER PATIENT MONITOR TRUSCOPE II



DESCRIPTION

Multi-Parameter Patient Monitor

This lightweight modular patient monitor is equipped with a 12.1" high resolution color TFT LCD display, offering the parameters ECG, RESP, NIBP, TEMP and SpO2.

FEATURES

- 12 lead ECG waveform display & total 11 waveforms display
- ST analysis, arrhythmia analysis
- Digital SpO2 even for low perfusion cases
- 72 hours tabular and graphical trends & data storage
- Drug dose and Hemodynamic calculation software