

MULTI PARAMETER PATIENT MONITOR TRUSCOPE CLASSIC



DESCRIPTION

15" Modular Multi-Para Patient Monitor

This modular patient monitor offers a wide range of standard and optional parameters.

FEATURES

- 15" high resolution color TFT LCD display with maximum of 10 waveforms display
- Basic 5 parameters - ECG, NIBP, SpO2, RR & Temperature
- ST analysis and arrhythmia analysis
- 96 hours tabular and graphical trends
- 7 lead ECG waveform display (12 lead optional)