

Suction Machine Apparatus YB-DX23B



OVERVIEW:

1. Adopt oil-impregnated vanes pump, high flow rate, efficient and durable.
2. One-way valve prevents positive pressure from entering into the work system.
3. Overflow protection device prevents liquid from entering into the inner pipe.
4. Stepless vacuum regulation is available based on the viscosity of sputum and mucus.
5. Hand switch and foot switch installed in parallel, can be selected at will.
6. Newly designed novel appearance with optimum height, easy to observe and exam from different position.

SPECIFICATION:

1. Max Vacuum: $\geq 0.09\text{MPa}$ (680mmHg)
2. Free airflow: $\geq 30\text{L/min}$
3. Liquid storage bottle: 2500ml \times 2(glass)
4. Power supply: AC220V 50Hz
5. Input: 280VA
6. Dimension: 465mm \times 425mm \times 920mm
7. G.W/N.W.: 24.5kg/22kg