

# MDX



## DESCRIPTION:

The GE Series of Anesthesia machines are a robust feat of engineering with heavy duty alloy frame and state-of-the-art design that ensures surgeons have one thing less to worry about during their work.

The compact design keeps the unit footprint to a minimal while providing all the features and technology found in most premium grade Anesthesia Machines today.

### **Special Feature**

Application: Adult, Pediatricλ

Driven mode: Pneumatically driven electronicallyλ  
controlled

Mode: Semi-open, Semi-closed or closed systemλ

Anesthetic vaporizer: One vaporizer (optional two)λ

### **Anesthesia Ventilator**

Display: High light LED displayλ

Ventilation Volume: 50~1600mlλ

Breath frequency: 1~90BPMλ

I:E ratio: 4:1 ~ 1:6λ

Monitor parameter: Tidal volume, breath frequency,  $MV\lambda$   
Alarm: MV high alarm, MV low alarm, TV high alarm, TV low alarm