



MAIN APPLICATIONS

Emergency Dept.

Home-Care

Patient Transfer

Room / Ward

Tracheostomy

Veterinary

Internal Rechargeable Battery

NEW ASKIR 230/12V BR is an electric suction unit for the aspiration of body liquids, oral, nasal and tracheal aspiration in adults or children. Suitable for emergency applications and for post-operative therapy, both in professional environment and home-care. Three different options for operation: AC/DC adapter, rechargeable battery and 12V make it a versatile suction unit approved for NON-STOP operation without overheating. Provided with visual and acoustic alarm indicating low battery level.

ACCESSORIES INCLUDED

- Liquid Collection Jar with overflow valve system (different options, see below)
- 6x10mm Silicone Tubes, length 140cm (autoclavable) + Conical Connector \varnothing 8-9-10mm
- Antibacterial & Hydrophobic Filter (single-patient) + CH20 Canula
- AC/DC Universal Adapter + Power Cord
- 12V Car Adapter

OPTIONS

- Stand on 5 Castors \varnothing 50mm with Brakes (antistatic upon request) REF 27730
- Carrying Bag (for 1000ml jar) SP 0207/02
- Carrying Bag (for 2000ml jar) SP 0207

TECHNICAL FEATURES

Motor	Oilless and maintenance-free piston pump
Power Feeding	14V \approx 4A by AC/DC Universal Adapter (model UE60-140429SPA1 (100-240V~ - 50 - 60Hz - 100VA)
	Internal Rechargeable Pb Battery (12V \approx 4A) - Autonomy 80min - Recharging Time 240min
	Car Adapter (12V \approx 1,9A)
ISO 10079-1 Classification	HIGH VACUUM / LOW FLOW
IP Code	IP21
Max Vacuum (adjustable)	-0.75 bar -75 kPa -563 mmHg
Max free air flow rate	16 l/min
Noise Level	63,0 dBA
Duty cycle	Non-stop operation (by mains only)
Weight	kg 3.50
Size	cm 35 x 21 x 18
Years of Warranty	2
Shipping carton	3
Place of Manufacturing	Italy

Italian
Medical
Touch
MADE IN ITALY

AVAILABLE CONFIGURATIONS

- RE 310211 NEW ASKIR 230/12V BR with 1000 ml autoclavable jar in polycarbonate Makrolon 2858
 - RE 310211/01 NEW ASKIR 230/12V BR with 2000 ml autoclavable jar in polycarbonate Makrolon 2858
 - RE 310211/03 NEW ASKIR 230/12V BR with 1000 ml Flovac ° disposable liner and reusable container
 - RE 310211/04 NEW ASKIR 230/12V BR with 2000 ml Flovac ° disposable liner and reusable container
- Jars in APEC ° autoclavable up to 143°C available upon request

CE 0123