



CLINEB is a professional piston compressor nebulizer for aerosoltherapy. Indicated for intensive treatments in clinics or hospitals, its design includes a large compartment for convenient storage of accessories and medications. Provided with pressure gauge and indicator.

STANDARD CONSUMABLES INCLUDED

- HI-FLO Jet Nebulizer
- Adult Masks with elastic strings
- Child Masks with elastic strings
- Mouthpiece
- Nosepiece (non invasive)

OPTIONS

- Air Tube
- Air Filter
- Detachable power cord

- Stand on Castors REF. 27731

TECHNICAL FEATURES

Compressor	Oiless and maintenance-free piston pump
Power Feeding	230V-50Hz - other configurations on request
Power consumption	184 VA
Max Pressure	2.5 bar 36 psi 250 kPa
Operating Pressure	1.30 bar 18.9 psi 130 kPa
Max Air Flow	16 l/min
Operating Air Flow	5,2 l/min
Neb-rate with 4 ml solution	0,50 ml/min
Mass Median Aerodynamic Diameter	3,25 μm (MMAD) according to EN-13544-1
Respirable Fraction (< 5 μm)	> 70%
Duty Cycle	Non-stop operation
Noise Level (EN 13544-1 Standard)	Approx. 57 dBA
Weight	2,20 kg
Size (Main Unit)	cm 25 x 19 x 16
Years of Warranty	2
Shipping carton	4
Place of Manufacturing	Italy

Italian
Medical
Touch
MADE IN ITALY

CODE

RE 300550 CLINEB

CE 0123