

Non-Contact Infrared Forehead Thermometer FS-300



Special Feature

- Quick Non-Contact Temperature Detection
- No Probe Filters, No Risk of Virus Transmission
- Versatility for General-Use Measurements
- Backlight function
- Automatically power off
- LCD display
- Pocket-size

Specifications

Measurement Method	Infrared measurement, Non-Contact measurement
Measurement distance	2 to 3cm from the measuring point
Measurement point	The measuring point is the temple area
Measurement range	Body temperature : 34°C to 42.5°C (93.2°F to 108.5°F) Accuracy ±0.2°C, (36°C to 39°C), Other Accuracy ±0.3°C Object temperature : 15°C to 60°C (59°F to 140°F) Accuracy ±2°C
Response time	Less than 2 seconds continuous measurement
Resolution	0.1°C
Memory	Automatically body temperature data (maximum 32) stored
Back-Light	If measurement complete, lighted 3 seconds
Automatic Power off	Auto Power-Off (after 1minute)
Accessory	FS-300 1set, Holder 1ea, Manual 1volume, Battery(AAA type) 2ea
Power supply	DC3V(AAA 1.5V Battery 2ea)